


INTERNATIONAL PRODUCT CATALOG



dentagine

SSWITE Dental®

Better Patient Outcomes
Improved Efficiency
Faster Practice Growth
















sswhitedental.com

Save 12 Hours Per Year in Chair Time Using the #1 Selling Metal Cutting Bur Worldwide

Innovative blade geometry combined with unparalleled neck strength produces a bur that quickly cuts metal with reduced chatter and breakage enabling one bur to remove up to three PFMs.

GREAT WHITE® GOLD SERIES















METAL CUTTING - CAVITY PREPARATION CARBIDE BURS

 <p>ROUND</p>	 <p>5 PK 50 PK</p> <p>GW2R 18202-5 300050</p> <p>GWSL2R 350146 10 PACK</p> <p>GW4R 18204-5 300051</p> <p>GWSL4R 350148 10 PACK</p> <p>GW6R 18206-5 300052</p> <p>GWSL6R 350150 10 PACK</p> <p>GW8R 18208-5</p>	 <p>PEAR</p>	 <p>5 PK 50 PK</p> <p>GW2P 18209-5</p> <p>GW245 18212-5 300053</p> <p>GW330 18210-5 300054</p> <p>GW330S 18211-5</p>
 <p>INVERTED CONE</p>	 <p>5 PK</p> <p>GW34 18214-5</p> <p>GW35 18215-5</p> <p>GW37 18216-5</p>	 <p>STRAIGHT FISSURE</p>	 <p>5 PK 50 PK</p> <p>GWH340 14918-5</p> <p>GWH342 14919-5</p> <p>GWH342L 14920-5</p> <p>GW1557 18221-5 300056</p> <p>GWSL1557 350167 10 PACK</p>
 <p>STRAIGHT FLAT END FISSURE</p>	 <p>5 PK 50 PK</p> <p>GW556 18217-5</p> <p>GW557 18218-5 300055</p> <p>GW557S 18219-5</p> <p>GWSL557 350162 10 PACK</p> <p>GW558 18222-5</p>	 <p>TAPER FLAT END FISSURE</p>	 <p>5 PK</p> <p>GW701 18223-5</p> <p>GW702 18224-5</p>
 <p>RESTORATION REMOVAL</p>	 <p>5 PK 50 PK</p> <p>GW1 (331) 18061-5 300057</p> <p>GW2 (1558) 18062-5 300059</p> <p>GW1S 18063-5 300058</p> <p>GW2S 18064-5 300060</p>	<p>TALON BUR</p>  <p>5 PK 50 PK</p> <p>T10 310297 300144</p>	<p><i>"I have yet to experience anything that cuts as fast, accurate and smooth as a Great White® Bur, saving my practice time and money while reducing stress. A single Great White® Bur replaces my use of two to three carbide burs for PFM crown removal, a time savings of about 12 hours per year and over \$4,000 in productivity."</i></p> <p><i>- Howard S. Glazer, DDS</i></p>


40% Faster Crown Preparations Versus Diamonds.

Great White Ultra Burs cuts both enamel and restorative materials faster than comparable diamond instruments, saving chair time and money.

GREAT WHITE® ULTRA CROWN & BRIDGE PREPARATION CARBIDE BURS

<p>FOOTBALL</p>  <p>10 PK HEAD LENGTH (mm) MAX HEAD DIAMETER (mm)</p>	 <p>GWU 379-023 350138 4.2 2.3</p>	<p>SHORT FLAT END TAPER</p>  <p>10 PK HEAD LENGTH (mm) MAX HEAD DIAMETER (mm)</p>	 <p>GWU 845KR-016 GWU 845KR-018 GWU 845KR-025 350140 350141 350142 4.0 4.0 4.0 1.6 1.8 2.5</p>	<p>DEPTH CUTTER</p>  <p>5 PK</p>	 <p>GWDC1 GWDC2 310158 310159</p>
<p>FLAT END TAPER</p>  <p>10 PK HEAD LENGTH (mm) MAX HEAD DIAMETER (mm)</p>	 <p>GWU 847-016 GWU 847-018 GWU 847-020 350144 350136 350137 8.0 8.0 8.0 1.6 1.8 2.0</p>	<p>FLAT END TAPER SHORT SHANK</p>  <p>10 PK HEAD LENGTH (mm) MAX HEAD DIAMETER (mm)</p>	 <p>GWU 847-016S GWU 847-018S GWU 847-020S 350131 350132 350130 8.0 8.0 8.0 1.6 1.8 2.0</p>		
<p>ROUND END TAPER</p>  <p>10 PK HEAD LENGTH (mm) MAX HEAD DIAMETER (mm)</p>	 <p>GWU 855-025 GWU 856-016 GWU 856-018 GWU 856-020 GWU856-023 350139 350143 350133 350134 350135 7.0 8.0 8.0 8.0 8.0 2.5 1.6 1.8 2.0 2.5</p>	<p>ROUND END TAPER SHORT SHANK</p>  <p>10 PK HEAD LENGTH (mm) MAX HEAD DIAMETER (mm)</p>	 <p>GWU 856-016S GWU 856-018S GWU 856-020S 350127 350128 350129 8.0 4.0 8.0 1.6 1.8 2.0</p>		

Crown and Bridge Preparation Kit



All-In-One Kit
Designed for the efficient atraumatic removal of failed amalgams, composites, zirconia and PFM Crowns.
Kit Reorder No. 18150

Kit Contents:
2 Great White® Ultra #856-016
2 Great White® Ultra #856-018
2 Great White® Ultra #856-020
1 Great White® Ultra #379-023
2 Great White® Ultra #847-016
2 Great White® Ultra #847-018
1 Autoclavable Bur Block

The Great White® Ultra:

- Cuts quickly and smoothly through enamel
- Negotiates amalgam and other restorative materials with minimal clogging and no "drag or stalling"
- Now a single instrument for completing crown preparations
- More efficient –less chair time
- Helps provide precise margins

"The Great White® Ultra is a versatile bur designed for crown and bridge. It cuts smoothly and efficiently through dentin, enamel, amalgam, composite and cast metal restorations. The unique blade geometry ensures smooth cutting action even when moving from dentin to amalgam or cast metal restorations." – Steven Abrams, DDS